



Gaumard®
Simulators for Health Care Education




HAL® S5301

Interdisciplinary and Immersive

“Hello, I’m HAL® — the world’s most advanced interdisciplinary patient simulator”

From emergency care to ICU and med-surg training, HAL is designed to fulfill educational objectives across clinical disciplines and blur the lines between simulation and real life. New Conversational Speech¹, lifelike motor movement, next-gen simulated physiology, UNI[®] 3, and many more industry-first capabilities usher in the next leap in simulation.



A new standard in patient simulation

HAL S5301 is the all-new flagship of the HAL brand, trusted by simulation programs worldwide for its innovation, ease of use, and reliability. The S5301 sets the standard for the future of patient simulation, building upon the design philosophy established with the original HAL launched nearly two decades ago. Join us as we chart HAL’s course into the future.





Introducing a new level of fidelity in neurological emergency simulation

HAL introduces the latest innovations in robotics and AI-enhanced speech¹ to simulate lifelike symptoms of a neurological emergency and enables team-based training through patient assessment and treatment without interruption. New features like facial droop, Conversational Speech¹, and arm motor control reproduce the progression of a stroke, helping teams train in time-dependent clinical skills and teamwork.

1. Conversational Speech¹ enhanced by artificial intelligence

HAL automatically answers medical history and initial assessment-related questions, follows voice commands, and learns over time¹.



2. Active motor function

Right arm motor reflex: shake hand, squeeze hand, raise arm, withdrawal response, and abnormal posturing. Head and eyes turn toward the provider speaking.

3. Abnormal eye movements

Consensual and nonconsensual pupillary response to light stimuli. Abnormal eye movements include strabismus, ptosis, and more.



4. Dynamic facial expressions

Lifelike facial droop and smiling, pained, quizzical, and scared facial expressions. Dynamic emotional states automatically express non-verbal cues including worry, anxiety, and lethargy.

5. Lifelike sensory response

Active pain response to pressure-sensitive sites: bilateral supraorbital notch, trapezius pinch (left shoulder), sternal rub, and right middle finger nailbed.

Hospital trauma team training made immersive

With new ultra-high fidelity anatomy and physiology, HAL supports advanced trauma care algorithms and essential surgical interventions using real tools and clinical techniques. Skin, bony landmarks, and internal tissue provide realistic tactile feedback to maximize the development of transferable skills. Internal sensors provide real-time feedback on provider interventions while automatically recording event data for debriefing.

1. Chest tube thoracostomy

Realistic left hemo/pneumothorax site supports palpation, incision, chest tube insertion, tube placement detection, bleeding, and suturing.

2. Surgical airway management

Anatomically accurate oral cavity and airway. Perform tracheotomy, cricothyrotomy, and retrograde intubation.

3. Abdominal bleeding wound

Penetrating abdominal wound responds to pressure and packing. Features internal, auto-refilling blood reservoir with 1.2-liter capacity.

4. Intraosseous access

Supports tibia and humeral intraosseous access and continuous infusion.

5. Real-time CPR performance feedback

Real-time quality feedback and reporting: Time to CPR, compression depth/rate, compression interruptions, ventilation rate, excessive ventilation, smart CPR coach.

6. Trauma arm and leg

Optional trauma arm and trauma leg accessories feature bleeding wound and tourniquet placement detection.





True-to-life imaging with Gaumard Ultrasound™

Facilitate ultrasound training through full-scale immersive emergency and trauma scenarios. Gaumard Ultrasound simulates the function and feel of a real portable ultrasound machine. Transducer range-of-motion is natural. Imaging is dynamic and lifelike, allowing students and professionals to refine imaging skills in team-based simulation.

The new HAL Emergency POCUS/eFAST module offers you a library of emergency ultrasound case imaging covering a wide variety of trauma presentations.



Learn more about
Gaumard Ultrasound at
www.gaumard.com/
HAL-S5301





The new gold standard in critical care in-situ simulation training

HAL S5301 revolutionizes critical care simulation training through powerful physiologic features and software. HAL helps participants improve skills and confidence by enabling learning experiences in real environments and using actual patient monitoring devices and mechanical ventilators.



Next-generation lung physiology and mechanical ventilator support

HAL S5301 features our latest advanced dynamic lung system design, capable of responding to mechanical ventilation with greater physiological accuracy. HAL interfaces with real mechanical ventilators like a real patient and supports standard ventilator modes, including AC, CMV, SIMV, and PSV, as well as PEEP and weaning strategies. The patented internal lung design means no manual calibration, proprietary adapters, or expensive and stationary external converters.

1. Real patient monitoring

Monitor vitals using real equipment: RR, 12-lead ECG, IBP, BP, SpO₂, TOF, and EtCO₂

2. High-quality auscultation

New, high-quality library of lung, heart, and bowel sounds. Auscultate using real stethoscopes.

3. Intra-arterial blood pressure monitoring

Radial arterial access site permits catheterization, flashback, sampling, and IBP monitoring; interfaces with real adjuncts, sensors, and devices.

4. Intravenous access

Features bilateral IV access sites, an antecubital vein blood draw site, and automatic virtual drug recognition on the lower left arm.

5. Blood glucose testing


Perform finger-stick glucose testing on the left index finger.

6. Train-of-Four monitoring

Supports Train-of-Four monitoring using real devices.

Urinary catheterization


Features internal 0.7-liter urine reservoir with variable urine/blood output.



Introducing the all-new UNI® 3

UNI 3 is our most powerful and intuitive patient simulator control software.

Manage vitals, track performance, and debrief with faster and easier-to-use tools designed to help you navigate even the most complex scenarios with ease.



Built new from the ground up


UNI 3 is built on a new, modern software platform, improving performance and stability with a refreshed yet familiar design.

New controls and tools

User account management stores personalized settings, scenario builder features visual guides to aid planning, and new diagnostics tools that ensure HAL is always operating optimally.

Feedback for debriefing

UNI 3 interfaces with sensors in HAL to provide you with more detailed real-time performance and event data to aid assessment and debriefing.




New Patient

First Name: Dorian
 Last Name: James
☐ In default patient

Info
 General Pain History Glasgow Coma Scale Labs

Street: 14070 SW 136 Street City: Miami
 Zip Code: 33196 State: Florida Country: USA
 Height: cm feet inch Weight: kg lbs
 Date of Birth: 5/15/1982 Birthplace: Miami
 Gender: Male Pronouns: he / him / o
 Relationship Status: Married Children:



Ready for training with the HAL® S5301 Simulation Learning Experiences™ package

The HAL S5301 Simulation Learning Experiences (SLE) package provides you with a library of 12 ready-to-use, evidence-based scenarios designed to help ease the facilitator's workload, increase realism, and standardize training.

The package includes accompanying UNI 3 preprogrammed scenarios that automatically manage the patient's vitals and a Facilitator's Guide book with comprehensive information for planning, setting up, and facilitating each SLE.

Myocardial infarction designer

Simulate a myocardial infarction with the easy-to-use 3D heart model and monitor the resulting changes using a real 12-lead ECG device or design your own rhythm using the point-by-point PQRST wave editor.

Patient profiles

The new Patient Chart lets you create simulated patients with detailed active and past medical histories. UNI® 3 seamlessly links medical chart data with Conversational Speech¹, enabling HAL to answer assessment and medical history questions automatically.

Scenario designer

Create scenarios tailored to your learning objectives and offer participants a wide range of standardized, repeatable learning events.

HAL® S5301 Select Features USA

General

- Height 5' 9" / 175 cm²
- Wireless and tetherless^{3,4}
- Microsoft Surface Pro preloaded with UNI® 3 Unified Simulator Control Software
- HAL S5301 Simulation Learning Experiences™ scenario package
- Bluetooth, Gaumard RF, and wired connectivity^{3,5}
- Compatible with Gaumard Ultrasound™
- Compatible with Care in Motion™
- Available in light, medium, or dark skin tone at no extra charge⁵

Neurological

- Active eye movement and object tracking
- Gaumard Conversational Speech¹
- Wireless streaming voice⁵
- Active neck movement and mouth movement
- Active facial expressions; left facial droop, right facial droop, pained, quizzical, scared, smiling
- Active emotional states: normal, worried, anxious
- Right arm motor reflex: shake hand, squeeze hand, raise arm, withdrawal response, and abnormal posturing (decorticate/ decerebrate)
- Stroke clinical presentations: facial droop, weakness in the right arm, abnormal posturing, and pain response
- Pressure sensitive sites: bilateral supraorbital notch, trapezius pinch (left shoulder), sternal rub, and right middle finger nailbed
- Train-of-Four monitoring using real devices
- Partial tonic-clonic and non-tonic-clonic seizures
- Sweating (diaphoresis) and tears

Airway

- Airway management: ETT, SGA, OPA, NPA⁷
- Tongue edema, laryngospasms, and pharyngeal swelling
- Tracheotomy, cricothyrotomy, and retrograde intubation

Breathing

- Spontaneous breathing
- Four anterior and posterior lung auscultation fields; new lung sound library
- Supports standard mechanical ventilators, modes of ventilation, weaning/liberation protocols; no calibration, proprietary adapters, or external converter adjuncts required
- Variable lung unit compliance, variable, bilateral and unilateral bronchi resistance, inspiratory effort and rate, respiratory drive, real CO₂ exhalation, Auto-PEEP
- Left hemo/pneumothorax
- Needle thoracentesis
- Presents capnography waveforms on real devices, including “shark fin” waveform

Cardiac

- Aortic, pulmonic, tricuspid, and mitral auscultation fields and new heart sound library
- Monitor vitals with real devices: 4-lead ECG, 12-lead ECG, pulse oximeter, NIBP/IBP monitor, capnograph
- ECG-derived respiratory monitoring
- UNI® 3D Myocardial Infarction Model
- UNI® 12-lead ECG waveform designer
- Defibrillation and pacing with live energy
- Anterior/lateral and anterior/posterior defibrillator pad placement
- Real-time CPR quality feedback
- Bilateral palpable pulses: carotid, brachial, radial, femoral, popliteal, pedal
- Circumoral skin coloration
- Bilateral IV access sites
- Radial arterial access site supports catheterization, flashback, sampling, and IBP monitoring
- Antecubital vein blood draw site on left arm
- Automatic virtual drug recognition at lower left arm
- Fingerstick glucose testing on the left index finger
- Tibia and humeral intraosseous access and infusion
- Capillary refill time testing at right middle finger

Genitourinary

- Male urinary catheterization with fluid return
- Computer-controlled urinary output rate and urine and/or blood mixture
- Internal, auto-refilling 0.7-liter urine reservoir

Gastrointestinal

- Four bowel auscultation fields and new bowel sound library
- Gastric distention

Trauma

- Internal, auto-refilling 1.2-liter blood reservoir
- Abdominal bleeding wound responds to pressure and packing
- Optional trauma arm and trauma leg accessories feature bleeding wound and tourniquet placement detection



Browse the full list of
HAL S5301 features at
[www.gaumard.com/
HAL-S5301](http://www.gaumard.com/HAL-S5301)

HAL® S5301 - Advanced interdisciplinary patient simulator

HAL® S5301 - Advanced Interdisciplinary Patient Simulator Package (USA)

S5301.PK

HAL S5301 patient simulator, Microsoft Surface Pro tablet preloaded with UNI 3.0, 4G LTE Wi-Fi router⁸, RF communications module, HAL S5301 Simulated Learning Experiences™ Scenario Package, Facilitator's Guide, abdominal wound insert, patient simulator accessories, two-day installation and training service, 1-Year Limited Warranty⁹. Extended service plans are available.

Gaumard Ultrasound™ System

30081159A

Package includes Gaumard Ultrasound laptop, transducer, Gaumard Ultrasound software license, and transport case.



HAL® S5301 Emergency Ultrasound POCUS/eFAST Pathologies Module

30081347A

HAL S5301 Emergency Ultrasound POCUS/eFAST Cases Module for Gaumard Ultrasound.

Traumatic Right Arm Amputation

30011856A

Lower-right arm with traumatic amputation and pressure-sensitive bleeding site.

Traumatic Left Leg Amputation

30011859A

Lower-left leg with traumatic amputation and pressure-sensitive bleeding site.

Gaumard Vitals™ Bedside Virtual Monitor

30081435A

Bedside, customizable virtual patient monitor. Package includes preconfigured all-in-one PC and one Gaumard Vitals patient simulator license.



Available in light, medium, or dark skin tone at no extra charge

Gaumard Vitals™ Portable Virtual Monitor

30081003A

Portable, customizable virtual patient monitor. Package includes preconfigured tablet PC and one Gaumard Vitals patient simulator license.

Care in Motion™ Mobile Video Debriefing System

CIM.PK

Care in Motion Tablet PC, 3 battery-powered HD wireless cameras, 3 adjustable camera grips, and transport case. One-Year Limited Warranty.

HAL[®] S5301

Interdisciplinary and Immersive



Gaumard[®]
Simulators for Health Care Education

Gaumard.com | Worldwide 305.971.3790 | Toll-Free USA & Canada 1.800.882.6655

1. Conversational Speech feature availability may vary by region or country. Internet access required. 2. Patient simulator approximate physical dimensions: height 5' 9" / 175 cm, weight 135lbs / 61kg. 3. Maximum wireless range will vary depending on environmental factors and conditions. 4. Battery is dependent on active features and settings; results may vary. 5. Some audio features are not available in long-range RF wireless mode. 6. Skin tone selection is available at the time of initial order only. 7. Dry exercises only; fluid insertion into the nasal and oral cavity is not supported. 8. Standard operation not dependent upon LTE plan. Contact your provider for details on LTE availability, support, compatibility, pricing, SIM card, and activation. 9. Warranty coverage, service, product installation, and training may not be available in all areas or countries. Product design and features are subject to change without notice. All trademarks and/or copyright materials are the property of their respective owners. Patented; other patents pending. © 2023 Gaumard Scientific, Inc. All Rights Reserved. 11190161A