



# ZOE<sup>®</sup> S504.200

## Gynecologic Skills Trainer

- Anatomically accurate pelvic landmarks
- Realistic normal and abnormal cervixes and uteri
- Modular additions for growing curricula
- Practice and perfect 'sensitive' procedures



## ZOE® Gynecologic Skills Trainer

S504.200.PK ● ● ●

The ZOE Gynecologic Torso comes with everything needed to start training in one complete package.

### Torso Features

- Adult-sized lower torso with ischial spines and sacrum
- Vaginal introitus facilitates placement of a female condom or diaphragm
- Interchangeable cervixes for visualization with speculum
- Lifelike cervixes and uteri support bimanual examination
- Palpable, realistic uteri with ovaries and fallopian tubes
- Minilaparotomy visualization and occlusion of fallopian tubes
- Realistic urethra and bladder for catheterization exercises
- Perform uterine sounding with real instruments
- Interchangeable vulva for optional ZOE packages

### Package Contents

- ZOE Gynecologic Torso
- Non-pregnant abdomen
- Anteverted uterus
- Retroverted uterus
- Clear IUD uterus
- Pregnant uteri: 6-8 weeks, 6-8 weeks w/ short ovarian ligaments, 10-12 weeks, 20 weeks
- (5) Normal patent cervixes
- (6) Abnormal cervixes
- Pregnant cervixes: (3) 6-8 weeks, (3) 10-12 weeks
- Urine kit
- Instruction manual

## ZOE® S504.200 Gynecologic Trainer

Lifelike anatomical features make the new ZOE the perfect addition to any gynecologic skills curricula, providing invaluable hands-on training experience on various gynecologic procedures.

- Speculum examination
- Bimanual examination
- Minilaparotomy
- Catheterization
- Fallopian occlusion
- Rectal medication
- Uterine sounding
- IUD placement

### Accurate anatomy

The new ZOE helps learners develop clinical skills for scenario-based exercises.

- Uniform skin with elasticity
- Ischial spines and coccyx
- Realistic cervixes & uteri

### Modular Additions

Optional additional pathology packages.

- Uteri external pathologies
- Uteri internal pathologies
- Fistulous tract models

## External Pathologies Uteri Package

S504.200.235

Interchangeable uteri for recognizing abnormalities.



### Package Contents

- Enlarged uterus
- Small uterus
- Uterus with moderate retroversion
- Myomatous uterus
- Uterus with left side salpingitis
- Uterus with right side salpingitis
- Severely anteverted-anteflexed uterus
- Uterus with large ovarian cyst
- Uterus with medium ovarian cyst
- Bicornuate uterus

## Internal Pathologies Uteri Package

S504.200.236

Interchangeable uteri for hysteroscopic viewing.



cutaway view

### Package Contents

- Normal anteverted uterus
- Uterus with polyposis
- Uterus with varied polyps
- Uterus with hyperplasia
- Myomatous uterus
- Uterus with early carcinoma
- Uterus with advanced carcinoma
- Uterus with fundus carcinoma
- Subseptate uterus

## Fistula Package

S504.200.700.1 ● ● ●

Lifelike practice suturing fistulous tracts.



### Package Contents

- Vulva with vesicovaginal, rectovaginal, and urethrovaginal fistulas
- Optional carrying case

## Request a Quote

[www.gaumard.com/quote](http://www.gaumard.com/quote)

[sales@gaumard.com](mailto:sales@gaumard.com)

Toll-Free USA & Canada

1.800.882.6655

Worldwide 305.971.3790