



Gaumard®
Simulators for Health Care Education




PEDIATRIC IV AND ARTERIAL ACCESS TRAINING ARMS

5-Year, 1-Year, and Newborn Skill Trainers

Advanced Pediatric IV and Arterial Access Training Arms.

The new Gaumard® Pediatric Training Arms offer educators the most realistic and versatile platform for teaching intravenous and arterial access skills to learners of all levels.



5-Year Pediatric IV and Arterial Access Training Arm

\$405.100 ● ● ●

Package includes: 5-Year Pediatric Training Arm, control box, blood dispensing bag, accessories, power supply, user manual, 1-Year Limited Warranty.

- Realistic arterial and venous network
- Adjustable pulse strength and rate (0 - 180 BPM)
- Palpable radial and brachial pulse
- Ulnar, cephalic, cubital, and basilic venipuncture
 - » infusion, bolus injection, blood draw/ collection, and observe flashback
- Radial and brachial arterial access
 - » blood sampling/ collection, and observe flashback
- IM injection
- Sub-Q injection
- Humeral IO access and continuous infusion
- External blood reservoir for continuous fluid supply and vessel pressurization
- Rotatable arm and shoulder
- Easily replaceable veins and arteries
- Lifelike skin
- Latex-Free
- Articulating wrist allowing palmar flexion
- Use as a stand-alone or attached to stand

1-Year Pediatric IV and Arterial Access Training Arm

\$406.100 ● ● ●

Package includes: 1-Year Pediatric Training Arm, control box, blood dispensing bag, accessories, power supply, user manual, 1-Year Limited Warranty.

- Realistic arterial and venous network
- Adjustable pulse strength and rate (0 - 180 BPM)
- Palpable radial and brachial pulse
- Ulnar and cubital venipuncture
 - » infusion, bolus injection, blood draw/ collection, and observe flashback
- Radial and brachial arterial access
 - » blood sampling/ collection, and observe flashback
- IM injection
- Humeral IO access and continuous infusion
- External blood reservoir for continuous fluid supply and vessel pressurization
- Rotatable arm and shoulder
- Easily replaceable veins and arteries
- Lifelike skin
- Latex-Free
- Articulating wrist allowing palmar flexion
- Use as a stand-alone or attached to stand

Newborn IV and Arterial Access Training Arm

\$408.100 ● ● ●

Package includes: Newborn Training Arm, control box, blood dispensing bag, accessories, power supply, user manual, 1-Year Limited Warranty.

- Realistic arterial and venous network
- Adjustable pulse strength and rate (0 - 180 BPM)
- Palpable radial and brachial pulse
- Ulnar and cubital venipuncture
 - » infusion, bolus injection, blood draw/ collection, and observe flashback
- Radial and brachial arterial access
 - » blood sampling/ collection, and observe flashback
- External blood reservoir for continuous fluid supply and vessel pressurization
- Easily replaceable veins and arteries
- Lifelike skin
- Latex-Free
- Articulating wrist allowing palmar flexion
- Use as a stand-alone or attached to stand
- Allows transillumination for vein access

Request a quote

www.gaumard.com/quote

sales@gaumard.com

Toll-Free USA & Canada

1.800.882.6655

Worldwide 305.971.3790